



Agence Nord:
ZA Object'ifs Sud - Lot A3
6 Allée Emilie du Châtelet
14123 Ifs
tél : 02.31.34.50.74
fax : 02.31.34.55.17



Agence Est:
Parc Club des Tanneries
2 Rue de la Faisanderie
67380 Lingolsheim
tél : 03.88.04.01.81
fax : 03.68.93.01.52

Agence Sud:
Bât Le Venango, 392 Rue Jean Dausset
AGROPARC - BP11575
84916 Avignon Cédex 9
tél : 04.90.27.17.95 fax : 04.90.27.17.52

www.deltalabo.fr
info@deltalabo.fr

Quick Thermal Conductivity Meter

QTM-710/700



SUMMARY

Quick and Easy Measurement of Thermal Conductivity of Sample!

The thermal conductivity of various samples such as block form or sheet in manufacturing, construction, handcraft, soil or the food industry can be measured with the Quick Thermal Conductivity Meter QTM-710/700.

Thermal conductivity can be measured quickly and easily and with good reproduction simply by placing the probe on the uniform-temperature sample surface.


- Thermal insulating material of fiber or foamed plastic
- Plastic, glass, wood, etc.
- Sheet, leather, etc.
- Bread dough, mashed food, powder, etc.
- High heat insulation material, ceramics, etc.



FEATURE

Clear and easy operation with the 5.7 inch color LCD touchscreen

The touchscreen enables intuitive operation. The history of the measured values can be called up by pressing the measurement display once, which quickly displays the necessary information and the operation menu.



Result window



Quick measurement with 3ch probes

Solid surfaces and film can be measured continuously without changing the probe when the same type or other types of various probe are connected.





Image of connection of various probes

Security function

Password settings and security level can be selected according to the intended use and environment. It is possible to restrict access to certain settings for non-management users.



Easy data transfer with USB

The measured data are easily transferred to personal computers using USB. The data are saved in CSV format and can be viewed and edited in spreadsheets.



The probe is traceability compliant

To further improve the reliability of the measured values a traceability system for the thermal conductivity of the standard substance (reference plate) has been set up.




FEATURE

Built-in thin film measurement

QTM-710

It is possible to measure the thermal conductivity of sheeted samples such as films or a sheet of fiber products.

Display of the average calculation of the measured value and temperature rise features enable easy analysis.



The sheeted sample can be measured without the use of a PC, and the detailed result can be printed out on a connected A4 printer.


Model:	QTM-710	Print: 2015/07/26 14:40
Serial:		Page: 1/1
Operator:	User	Date: 2015/11/19 12:19
Sample No:	15-001	Sample ID: TF 0.3mm PD-11 6 Q3
Meas.Mode:	Thin Film Measurement	
Probe Unit:	PD-11N	
Probe No:	00006	
Probe Const:	Main	
λ :	0.2035 [W/(m K)]	
	λ [W/(m K)]	1st 2nd 3rd 4th 5th AVG.
Polye No.6 QTM70003	0.03637	0.8303 0.8394 0.8342 0.8294 0.8451 0.8357
SiliconNo.6 QTM70003	0.2228	-0.03954 -0.02269 -0.03809 -0.01906 -0.03984 -0.03184
QuartzNo6 QTM70003	1.425	-0.3844 -0.4255 -0.4093 -0.3787 -0.4767 -0.4149

Box type probe PD-11N

QTM-710

QTM-700

Thermal conductivity is measured by simply placing the probe on a smooth sample surface. The PD-11N enables various measurements of the sample in block, powder and sheet.



- Heat insulation material
- Ceramics
- Rubber
- Plastic
- Glass

Insulated moisture-proof probe PD-13N

QTM-710

QTM-700

PD-13N is covered the insulation film to PD-11N. Thermal conductivity of aqueous materials and conductive materials are possible to measure with PD-13N.



- Food
- Concrete paste
- Metal

High-temperature probe PD-31N

QTM-710

QTM-700

Thermal conductivity is measured by placing the sample between the heater and the temperature sensor. The PD-31N is a conventional type of sensor and allows the thermal conductivity in higher temperature ranges to be measured.



- New material
- Firebrick
- Heat insulation material


Vessel for powder

QTM-710

QTM-700

Thermal conductivity in the bulk density *of the powder sample can be measured using the optional container for powder (using scales). (The mass of the collected sample must be known.)

*To calculate the bulk density the mass of the sample is required.




MEASUREMENT PRINCIPLE

Hot Wire Method (Transient Hot Wire Method)

The temperature of the heater wire will rise since the heat is generated when applying the constant current to the heater wire stretched in a straight line inside the sample. The temperature rise will be a straight line as shown in the right figure when the time axis is switched to a logarithmic scale.

The angle of this straight line is in reverse proportion to the thermal conductivity of the sample. The angle becomes large if the thermal conductivity of the sample is low, and the angle becomes small if it is high. Therefore, the thermal conductivity can be obtained from this angle.

QTM-710/700 measures the thermal conductivity with the box type probe (PD-11N and PD-13N) which uses the hot wire method and with the high-temperature probe (PD-31N) based on the hot wire method.



Thin film measurement

The measurement of thin sheet samples for thermal conductivity is based on the comparison of the sample with a reference plate during a temperature rise when both are heated by the probe.

Before measuring, prepare two to four homogeneous samples of known thermal conductivity. (Fig.1.)




Fig. 1 Measuring method


Place a sample closely on the plate and heat it with the probe. Then, monitor the temperature rise versus logarithmic time for some period to find the deviation from the initial state. Apply this method to a couple of plates and obtain the relation between the thermal conductivity and the deviation.

When the thermal conductivity of the sample is higher than that of reference plate, the deviation becomes positive value, and when it is lower, the deviation becomes negative.

Find the point of intersection where the deviation equals zero on the ε (deviation) versus λ (thermal conductivity) plot.

Measurement example (thin film measurement)

PTFE sheet



Ref	Ref λ	1	2	3	4	5	AVG
Polye No.6 QTM	0.03637	0.8303	0.8394	0.8342	0.8294	0.8451	0.8357
Silicone No.6 QTM	0.2228	-0.03954	-0.02269	-0.03809	-0.01906	-0.03984	-0.03184
Quartz No.6 QTM	1.425	-0.3844	-0.4255	-0.4093	-0.3787	-0.4767	-0.4149

PRODUCT STRUCTURE / SPECIFICATION

Product structure

■ QTM-710 <for sheet sample>

- Main unit
- Aluminum cooling plate
- Reference plate (clear quartz/silicone rubber/polyethylene foam)
- Probe is selectable (PD-11N or PD-13N)

■ QTM-700 <Not for sheet sample>*

- Main unit
 - Aluminum cooling plate
 - Reference plate (silicone rubber)
 - PD-11N or PD-13N
- ※Note that the thin sheet sample is unable to measure with QTM-700.

Specification

Specification	Contents	
Name	Quick Thermal Conductivity Meter	
Type	QTM-710	QTM-700
Method	Hot Wire Method	
Accuracy※ ¹	Within $\pm 5\%$ reading value per reference plate※ ² (room temp. only)	
Repeatability※ ¹	3% (at the measurement of reference plate)	
Thin film measurement	Able	Unable
On-window display	5.7 inch color LCD	
Language	English / Japanese / Mandarin Chinese / Korean	
External I/O	RS232C : 2ch (Printer, Data Acquisition Software) USB : 2ch (Printer, USB flash drive)	
Ambient conditions	Temperature : 5 to 35°C Humidity : 85%RH or below (no condensation)	
Power supply	DC24V 5A (main unit), AC100-240V $\pm 10\%$ 50/60Hz (AC adaptor)	
Dimensions	262 (W) \times 276 (D) \times 158 (H) mm	
Weight	Approx. 4kg	
CE marking	CE marking (EN61326-1, EN61010-1), suitable for RoHS	

※¹ Accuracy and repeatability of thin film measurement depend on the state of the sample.

※² Reference plate (clear quartz / silicone rubber / polyethylene foam)

<Normal measurement>

Probe	Box type probe PD-11N	Insulated moisture-proof probe PD-13N	High-temperature probe PD-31N
Method	Hot Wire Probe Method		Hot Wire Method
Measurement range	0.03~12W/(m K)		0.06~5W/(m K)
Dimensions	110 (W) \times 50 (D) \times 100 (H) mm		
Measurement temperature	5~35°C		-100~1000°C
Measurement time	60sec		60sec※
Sample size	100 \times 50 \times 20 mm or more		100 \times 50 \times 20 mm or more/2 pieces

※Optionally set the measurement time.

<Thin film measurement>

Probe	Box type probe PD-11N	Insulated moisture-proof probe PD-13N
Measurement range	0.03~5W/(m K)	
Measurement temperature	5~35°C	
Measurement time	60sec	
Sample size	Size: 100 \times 50 mm or more Thickness: Indication as to the thermal conductivity of the sample. 0.1W/(m K) \rightarrow 0.1mm or more 5 W/(m K) \rightarrow 5mm or more	

Part Number	Part Name	Remarks
12-06477	PD-11N Box Type Probe	
12-06477-01	PD-13N Insulated Moisture-proof Probe	
12-06479	PD-31N High-temperature Probe	
12-06753	Reference Plate/Clear Quartz for QTM-700	Approx. $\lambda = 1.4 \text{ W/(m K)}$
12-06753-01	Reference Plate/Silicone Rubber for QTM-700	Approx. $\lambda = 0.2 \text{ W/(m K)}$
12-06753-02	Reference Plate/Polyethylene Foam for QTM-700	Approx. $\lambda = 0.03 \text{ W/(m K)}$
12-06753-03	Plate/Silicone Sponge H for QTM-700	Approx. $\lambda = 0.07 \text{ W/(m K)}$
12-06753-04	Plate/Silicone Sponge L for QTM-700	Approx. $\lambda = 0.1 \text{ W/(m K)}$
12-06753-05	Plate/Zirconia for QTM-700	Approx. $\lambda = 3 \text{ W/(m K)}$
12-06753-06	Plate/Mullite for QTM-700	Approx. $\lambda = 5 \text{ W/(m K)}$
12-06753-07	Plate/Hastelloy C for QTM-700	Approx. $\lambda = 9 \text{ W/(m K)}$
20-07341	Aluminum Cooling Plate	
12-01876	Vessel for Powder	for PD-11N
69-00670	Brush	
IDP-100	Dot Matrix Printer	
64-00625	Connection Cable (MiniDIN8P-DSUB9PM) 160mm	for IDP-100/SOFT-CAP
12-02013	RS-232C Connection Cable (9P-25P) 1m	for IDP-100
12-04260	Printing Roll RP5860 4rolls Set	for IDP-100
69-00719	Ribbon Cartridge IR-91B Black	for IDP-100
DP-600	Thermal Printer	
64-00643-03	USB Cable	for DP-600
69-00522-01	Thermal Roll Paper P-58-30 (10 rolls)	for DP-600
12-03265	Data Acquisition Software (SOFT CAP)	
12-02012	RS-232C Connection Cable (9P-9P) 2m	for SOFT-CAP
64-00177-00	USB Serial Transducer US232R-10	for SOFT-CAP


The above reference plates and plates are also common to QTM-710.



Thermal Printer
DP-600



Dot Matrix Printer
IDP-100



Vessel for powder
12-01876



Reference plate/Aluminum cooling plate

KEM KYOTO ELECTRONICS
MANUFACTURING CO.,LTD.
<http://www.kyoto-kem.com>

Overseas Division : 2-7-1, Ichigaya-sadohara-cho, Shinjuku-ku
TOKYO, 162-0842, JAPAN
Fax : +81-3-3268-5591 Phone : +81-3-5227-3156

Distributed by

Agence Nord:

ZA Object'ifs Sud - Lot A3
6 Allée Emilie du Châtelet
14123 Ifs
tél : 02.31.34.50.74
fax : 02.31.34.55.17

Agence Sud:

Bât Le Venango. 392 Rue Jean Dausset
AGROPARC - BP11575
84916 Avignon Cédex 9
tél : 04.90.27.17.95 fax : 04.90.27.17.52



Agence Est:

Parc Club des Tanneries
2 Rue de la Faisanderie
67380 Lingolsheim
tél : 03.88.04.01.81
fax : 03.68.93.01.52

www.deltalabo.fr
info@deltalabo.fr