

Kinematic Viscosity of Asphalts (Bitumen)

ASTM D2170 - EN/ISO 12595



Completely stainless steel
Conforms to ASTM D2170 & D445
High precision stability $\pm 0.01^\circ\text{C}$
Large windows
Easy to operate
RS232 communication

General

ASTM D2170 method covers procedures for the determination of kinematic viscosity of liquid asphalts (bitumens), road oils and distillation residues of liquid asphalts (bitumens) all at 60°C (140°F) and of asphalt cements at 135°C (275°F).

Construction

The Tamson TV2000 and TV4000 viscometer baths are specially designed for tests that require ultra-precise temperature control, or processes that need to be followed visually. The stainless steel construction with 25 mm thick rock wool insulation ensures exceptionally stable temperatures which is further improved by an ingenious stirring mechanism with baffle plates. All wetted parts are made of stainless steel, providing resistance against all usual bath fluids. The bath is fitted with adjustable feet for leveling. The cover of the bath has a number of round $\varnothing 51$ mm openings with lids, for suspending glass capillary viscometers in holders. To work at temperatures lower than ambient $+5^\circ\text{C}$, use of cooling must be made.

Item	Unit	TV2000	TV4000
P/N 230V/50-60Hz		00T0782	00T0772
P/N 115V/60Hz		00T0784	00T0774
Range	Ambient..230°C/446°F		
Reading	°C or °F menu selectable		
Interface	RS232		
Setting ±	[°C]	0.01	
Stability* ±	[°C]	0.01	
Uniformity* ±	[°C]	0.01	
Heating 230V	[kW]	2.8	2.8
Heating 115V	[kW]	2.3	2.3
Heaters		2	2
Bath volume	[L]	20	40
Number of lids		3 or 4	7 or 8
Opening lid	[mm]	ø 51	ø 51
Window	[mm]	140 * 285	270 * 285
Opening bath	[mm]	130 * 165	260 * 240
Depth	[mm]	300	300
Length	[mm]	350	400
Width	[mm]	470	590
Height	[mm]	590	590
Weight	[kg]	40	41
Power	[W]	3100	3100
Frequency	[Hz]	Suited for both 50 & 60	
CE	Conforms to CE regulation		
* Value measured in water @ +50°C/ +122°F			

Cooling fluid can be pumped through the cooling coil inside the apparatus. Tap water or a combination with the external Tamson TLC10-3 or TLC15-5 cooling circulator can be used for this purpose. All models are fitted with double windows in front and rear walls. The windows are formed with two panes of tempered safety glass separated by 20 mm of air space. Visibility through the bath is excellent. A power plug on the backside is mounted to provide power for an optional Z41 LED illumination unit.

Agitation

A vane type stirrer with brass bearings move the bath fluid past the heaters and then from under the main baffle plate, thus directing the freshly heated fluid inside bath.

Span

All baths can be operated from ambient $+5$ up to $+230^\circ\text{C}$ ($41..446^\circ\text{F}$). With the use of the built-in cooling coil, span lies 5°K above the temperature of the cooling liquid.

Kinematic Viscosity of Asphalts (Bitumen)

ASTM D2170 - EN/ISO 12595

Accuracy and set point

The set point can be set in steps of 0.01°. The system overall accuracy is within $\pm 0.01^{\circ}\text{C}$. After the temperature control is stable, the offset can even be adjusted with $\pm 0.005^{\circ}\text{C}$.

Viscometer arrangement

A cover with seven openings with lids, arranged in two rows of resp. four and three, is provided for the TV4000. The TV2000 offers three openings (2+1). These $\varnothing 51$ mm openings will accommodate glass capillary viscometer in holders (please see our next page). Additionally two separate thermometers can be placed through two $\varnothing 12.5$ mm openings in the cover.

Safety

The bath conforms to CE regulations. It also is equipped with a mechanical adjustable and resettable safety thermostat. Advanced safety features are microprocessor control of:

- Electronic- and processor system
- Control and feedback from each heating element
- System accuracy

System error results in total cut off from the power supply.

Optional equipment specifically for ASTM D2170
Please see the next page.

TV2000 is standard included with:		
P/N	Picture	Description
23T2404		Cover with 3 openings: - 3 x $\varnothing 51$ mm opening - 2 x $\varnothing 12.5$ mm opening for thermometer
		3 * lid for $\varnothing 51$ mm opening

Optional cover for TV2000:		
P/N	Picture	Description
23T2405		Cover with 4 openings: - 4 x $\varnothing 51$ mm opening - 2 x $\varnothing 12.5$ mm opening for thermometer
		4 * lid for $\varnothing 51$ mm opening

TV4000 is standard included with:		
P/N	Picture	Description
23T2400		Cover with 7 openings: - 7 x $\varnothing 51$ mm opening - 2 x $\varnothing 12.5$ mm opening for thermometer
		7 * lid for $\varnothing 51$ mm opening

Kinematic Viscosity of Asphalts (Bitumen)

ASTM D2170 - EN/ISO 12595

Optional cover for TV4000:		
P/N	Picture	Description
23T2401		Cover with 8 openings: - 8 x ø51 mm opening - 2 x ø12.5mm opening for thermometer
		8 * lid for ø 51 mm opening

Accessories for ASTM D2170:			
P/N	Picture	TV2000	TV4000
00T0908		Illuminator "Z41" backpanel (85~230V/50-60Hz)	
10T6090		Timer, 8 positions	
E20 Thermometer		Please see specification sheet "E20 thermometers"	
14T0303		Adapter to insert an E20 thermometer in the opening of the cover	
Cannon Fenske Opaque (CFO) viscometer, Zeitfuchs Cross Arm® viscometer (ZCA), BS/IP/RF viscometer		Please see specification sheet "Viscosity accessories".	
06T1724		Stopper for viscometer CFO/ZCA (pack of 12 pieces)	
25T1069		Special stopper for BS/IP/RF viscometer	
10T6071		Viscometer holder for Cannon Fenske Opaque	
10T6327		Viscometer holder for Zeitfuchs Cross Arm®	

Tamson Instruments

Specification sheet

Kinematic Viscosity of Asphalts (Bitumen)

ASTM D2170 - EN/ISO 12595

Accessories needed for ASTM D2170:			
P/N	Picture	TV2000	TV4000
10T6051		Viscometer holder for BS/IP/RF	
08T0001		Silicon oil, suitable from 20 to 150°C (1 can of 20L needed)	Silicon oil, suitable from 20 to 150°C (2 cans of 20L needed)
Reference standards		Please see specification sheet "Viscosity accessories".	

Bath options:			
P/N	Picture	TV2000	TV4000
07T0086		Level detector/float	
07T0087			Level detector/float
02T0203		Spill tray	
02T0201			Spill tray

Agence Nord:

Bâtiment LETNA
Boulevard de l'espérance
14123 Cormelles le Royal
tél : 02.31.34.50.74
fax : 02.31.34.55.17



Agence Sud:

Bât Le Venango. 392 Rue Jean Dausset
AGROPARC - BP11575
84916 Avignon Cédex 9
tél : 04.90.27.17.95 fax : 04.90.27.17.52

Agence Est:

Parc Club des Tanneries
2 Rue de la Faisanderie
67380 Lingolsheim
tél : 03.88.04.01.81
fax : 03.68.93.01.52

www.deltalabo.fr
info@deltalabo.fr