

Agence Nord:
ZA Object'Ifs Sud - Lot A3
6 Allée Emilie du Châtelet
14123 Ifs
tél : 02.31.34.50.74
fax : 02.31.34.55.17



Agence Est:
Parc Club des Tanneries
2 Rue de la Faisanderie
67380 Lingolsheim
tél : 03.88.04.01.81
fax : 03.68.93.01.52

Agence Sud:
Bât Le Venango. 392 Rue Jean Dausset
AGROPARC - BP11575
84916 Avignon Cédex 9
tél : 04.90.27.17.95 fax : 04.90.27.17.52

www.deltalabo.fr
info@deltalabo.fr

51-21... Geltimer GT-S

Standards

DIN 16916, DIN 16945, DIN EN ISO 9396



Application

The Geltimer GT-S is intended for the measurement of Gelation time of reaction resins in accordance with DIN 16 945 Page 1, DIN 16 919 and DIN EN ISO 9396.

Features

Measurement on the Geltimer GT-S is performed under controlled conditions, is automatic and extremely simple to perform. It features exceptionally good reproducibility. Measurements can be made on any desired reaction resin. The instrument is extremely easy to operate and functions virtually maintenance-free.

A measuring stamper of aluminium, or steel for aggressive materials, is moved down and up in the reaction resin which is held in a test tube. When the reaction point is reached the test tube with the resin mixture and the measuring stamper is drawn upward in the reciprocating motion and an optical sensor stops the clock which was set in motion at the commencement of the test. The gelation time can then be read off directly on the clock.

As a result of the modular layout the lower part of the housing can be extended to incorporate electric heating or a thermostatically controlled oil or water bath.

Technical Data (Baseunit 51-210-014)

| | |
|-----------------------|--|
| Time Measurement | 0 to 99 h 59 min 59 s |
| Accuracy | ± 1 second / 2 hours |
| Reciprocating cycle | every 10 secs. (On request other times can be offered) |
| Operating temperature | 0 to 50 °C |
| Storage temperature | -30 to 70 °C |

Technical Data (Heating Block 51-210-018 + 51-210-019)

| | |
|----------------------|--|
| Temperature range | 0 to 250 °C (Actual value display) 30 to 200 °C (on request 250 °C) – (Set value) |
| Temperature accuracy | ± 0,3 % of the scale range ± 1 digit |
| Resolution | 0,1 °C |
| Stability | ± 0,2 °C |



Dimensions and Connection (Baseunit 51-210-014)

Dimensions (ØxH) 130 x 420 mm
 Weight ~ 2,8 kg
 Mains 230 VAC / 24 VAC, 50 / 60 Hz
 Power 5 VA

Dimensions and Connection (Heating Block 51-210-018 + 51-210-019)

Dimensions (Ø) 130 mm
 Mains 100 VAC ... 240 VAC, 50 / 60 Hz
 Power 500 Watt

Partlist

| Item no. | Description |
|---------------|---|
| 51-210-014 | Geltimer GT-S |
| 51-210-018 | Heating Block Ø 16mm |
| 51-210-019 | Temperature controller for heating block |
| 51-210-013 | PE foam ring (10 pieces) |
| 51-216-001 | Stampers made of aluminium for Geltimer (500 Stück) |
| 0400-22-16160 | Test tubes (100) 16 x 160 mm |

Agence Nord:

ZA Object'ifs Sud - Lot A3
 6 Allée Emilie du Châtelet
 14123 Ifs
 tél : 02.31.34.50.74
 fax : 02.31.34.55.17

Agence Sud:

Bât Le Venango. 392 Rue Jean Dausset
 AGROPARC - BP11575
 84916 Avignon Cédex 9
 tél : 04.90.27.17.95 fax : 04.90.27.17.52



Agence Est:

Parc Club des Tanneries
 2 Rue de la Faisanderie
 67380 Lingolsheim
 tél : 03.88.04.01.81
 fax : 03.68.93.01.52

www.deltalabo.fr
info@deltalabo.fr