



Agence Nord :
Plug N'Work - Campus Effiscience
2 Rue Jean Perrin - Bât D
14460 Colombelles
Tél: 02.31.34.50.74 Fax: 02.31.34.55.17

Agence Sud :
Créativa Bât A - Site Agroparc
200 Rue Michel de Montaigne - 84911 Avignon
Tél : 04.90.27.17.95 Fax : 04.90.27.17.52

www.deltalabo.fr

SETA COMPACT CLOUD & POUR POINT CRYOSTAT

Convenient refrigeration of fuel samples in accordance with all primary cloud and pour point test methods.

ASTM D97-IP 15; ASTM D2500-IP 219; ASTM D5853 - IP441; ISO 3015; ISO 3016; EN 23015; NF T60-105; FTMS 791b 201; JIS K2269



- TEMPERATURE RANGE 0°C TO -34°C
- 12 TEST POSITIONS
- ENERGY EFFICIENT OPERATION
- SPACE SAVING, SMALL FOOTPRINT

The Seta 'Compact' Cryostat has a 2 stage CFC free refrigeration system. The unit has 3 individually controlled test compartments, each can accommodate up to 4 sample wells.

The bath is supplied with 12 test jars (each comprising graduated glassware, thermometer cork, insulating gasket and disk), a compartment thermometer, and instruction manual.

Seta Product Ref	94100-0
Temp range:	0°C to -34°C
Ambient temp range requirement:	15° to 32°C
Cooldown time:	Approx 1hr from 32°C to -34°C
Compartments	3 individually controlled
Test positions	12
Bath liquid capacity:	Approx 2 litres
Refrigerant Gases:	CFC free First stage R134a; Second stage R507
Integral recharge:	'Spare' R507 canister built in
Power:	750W
Voltage:	110/120V 60Hz or 220/240V 50Hz
Size (HxWxD):	60 x 60 x 85cm
Weight:	115kg