

## Fuels Testing Standards

Product Code	Test Method	Product Nature	Certified Value	Pack Size
CRM-ADKR	Acidity ASTM D 3242	Kerosene	0.0021 mg KOH/g	250mL
CRM-APKR	Aniline Point ASTM D 611	Jet Aviation fuel	62.30°C	250mL
CRM-ACKR	Aromatics Content ASTM D 5580	Gasoline	28.25%	250mL
CRM-BEGA	Benzene Content ASTM D 4053 / D 5580	Gasoline	0.827 % Vol.	250mL
CRM-CNGO	Cetane Number ASTM D 613	Diesel	52.6	1000mL
CRM-FIKR	FIA Aromatics ASTM D 1319	Kerosene	15.18 % vol	250mL
CRM-OMGA	Motor Octane Number ASTM D 2700	Gasoline (Unleaded)	86.9	1000mL
CRM-ORGA	Res. Octane Number ASTM D 2699	Gasoline (Unleaded)	98.0	1000mL
CRM-VPGA	Reid Vapour Pres. ASTM D 4953A, 323A	Gasoline	82.6 kPa	250mL

*All values stated are nominal - 12/07/2011*