

KDS101

Syringe	Diameter	Minimum	Maximum
10 µl	0.46 mm	0.001 µl/min	0.350 µl/min
25 µl	0.73 mm	0.001 µl/min	0.884 µl/min
50 µl	1.03 mm	0.001 µl/min	1.759 µl/min
100 µl	1.46 mm	0.001 µl/min	3.526 µl/min
250 µl	2.3 mm	0.01 µl/min	8.78 µl/min
500 µl	3.26 mm	0.01 µl/min	17.65 µl/min
1 ml	4.61 mm	0.1 µl/min	35.2 µl/min
2.5 ml	7.28 mm	0.1 µl/min	88.0 µl/min
3 ml	8.59 mm	0.1 µl/min	122.5 µl/min
5 ml	10.3 mm	0.1 µl/min	176.2 µl/min
10 ml	14.57 mm	0.001 µl/min	0.351 ml/min

	DELTA LABO Créativa Bât A - Site Agroparc 200 Rue Michel de Montaigne BP 21221 84911 AVIGNON Cédex 9
VOTRE PARTENAIRE LABORATOIRE	Tél / fax 04 90 22 04 27 www.deltalabo.fr