

KDS120

| Syringe | Diameter | Minimum | Maximum |
|----------------|-----------------|----------------|----------------|
| 10 µl | 0.46 mm | 0.1 µl/hr | 126.5 µl/hr |
| 25 µl | 0.73 mm | 0.1 µl/hr | 318.8 µl/hr |
| 50 µl | 1.03 mm | 0.2 µl/hr | 625 µl/hr |
| 100 µl | 1.46 mm | 1 µl/hr | 1274 µl/hr |
| 250 µl | 2.3 mm | 2 µl/hr | 3164 µl/hr |
| 500 µl | 3.26 mm | 3 µl/hr | 6359 µl/hr |
| 1 ml | 4.61 mm | 0.01 ml/hr | 13.2 ml/hr |
| 2.5 ml | 7.28 mm | 0.02 ml/hr | 31.7 ml/hr |
| 3 ml | 8.59 mm | 0.02 ml/hr | 44.9 ml/hr |
| 5 ml | 10.3 mm | 0.03 ml/hr | 87.0 ml/hr |
| 10 ml | 14.57 mm | 0.1 ml/hr | 125.0 ml/hr |

| | |
|---|---|
|  | DELTA LABO Créativa Bât A - Site Agroparc 200 Rue Michel de Montaigne BP 21221 84911 AVIGNON Cédex 9 |
| VOTRE PARTENAIRE LABORATOIRE | Tél / fax 04 90 22 04 27 www.deltalabo.fr |