

KDS200/220, KDS210/230, KDS250/260, KDS270, KDS330

| Syringe | Diameter | Minimum | Maximum |
|----------------|-----------------|----------------|----------------|
| 10 µl | 0.46 mm | 0.001 µl/hr | 21.1 µl/min |
| 25 µl | 0.73 mm | 0.003 µl/hr | 53.15 µl/min |
| 50 µl | 1.03 mm | 0.005 µl/hr | 105.8 µl/min |
| 100 µl | 1.46 mm | 0.009 µl/hr | 212.6 µl/min |
| 250 µl | 2.3 mm | 0.021 µl/hr | 527.6 µl/min |
| 500 µl | 3.26 mm | 0.042 µl/hr | 1060 µl/min |
| 1 ml | 4.61 mm | 0.083 µl/hr | 2119 µl/min |
| 2.5 ml | 7.28 mm | 0.207 µl/hr | 5286 µl/min |
| 3 ml | 8.59 mm | 0.288 µl/hr | 7360 µl/min |
| 5 ml | 10.3 mm | 0.414 µl/hr | 634 ml/hr |
| 10 ml | 14.57 mm | 0.828 µl/hr | 1270 ml/hr |
| 20 ml | 19.05 mm | 1.414 µl/hr | 2171 ml/hr |
| 30 ml | 21.59 mm | 1.817 µl/hr | 2789 ml/hr |
| 50 ml | 28.9 mm | 3.277 µl/hr | 4998 ml/hr |
| 60 ml | 26.6 mm | 2.757 µl/hr | 4234 ml/hr |
| 100 ml | 34.9 mm | 4.746 µl/hr | 7289 ml/hr |
| 140 ml | 38.4 mm | 5.746 µl/hr | 8824 ml/hr |

