

**KDS270**

<b>Syringe</b>	<b>Diameter</b>	<b>Minimum</b>	<b>Maximum</b>
10 µl	0.46 mm	0.001 µl/hr	21.1 µl/min
25 µl	0.73 mm	0.003 µl/hr	53.15 µl/min
50 µl	1.03 mm	0.005 µl/hr	105.8 µl/min
100 µl	1.46 mm	0.009 µl/hr	212.6 µl/min
250 µl	2.3 mm	0.021 µl/hr	527.6 µl/min
500 µl	3.26 mm	0.042 µl/hr	1060 µl/min
1 ml	4.61 mm	0.083 µl/hr	2119 µl/min
2.5 ml	7.28 mm	0.207 µl/hr	5286 µl/min
3 ml	8.59 mm	0.288 µl/hr	7360 µl/min
5 ml	10.3 mm	0.414 µl/hr	634 ml/hr
10 ml	14.57 mm	0.828 µl/hr	1270 ml/hr
20 ml	19.05 mm	1.414 µl/hr	2171 ml/hr
30 ml	21.59 mm	1.817 µl/hr	2789 ml/hr
50 ml	28.9 mm	3.277 µl/hr	4998 ml/hr
60 ml	26.6 mm	2.757 µl/hr	4234 ml/hr

	<b>DELTA LABO</b> Créativa Bât A - Site Agroparc 200 Rue Michel de Montaigne BP 21221 84911 AVIGNON Cédex 9
	VOTRE PARTENAIRE LABORATOIRE Tél / fax 04 90 22 04 27 <a href="http://www.deltalabo.fr">www.deltalabo.fr</a>