


## High Temperature Viscosity Standards



High Temperature Type S3S

Overview	Specification	Shipping	Related Products
<p>Paragon ISO 17025 / ISO Guide 34 certified high temperature viscosity standards are for the calibration and verification of glass capillary viscometers and other viscosity measuring equipment where operating temperature is controlled precisely.</p>			
 No 0649		 No 4589	

Products	Kinematic mm <sup>2</sup> /s (cSt)		Dynamic mPa.s (cP)			Density g/mL		SUS	SFS		
Product Code	20°C 68°F	25°C 77°F	37.78°C 100°F	40°C 104°F	50°C 122°F	60°C 140°F	80°C 176°F	98.89°C 210°F	100°C 212°F	135°C 275°F	150°C 302°F
S3S	4.596	4.059	3.058	2.926	2.423	2.047	1.530	1.212	1.197	-	0.7424
S6S	10.59	8.945	6.113	5.761	4.500	3.618	2.502	1.880	1.852	-	1.066
S20S	43.32	33.89	19.63	18.02	12.70	9.360	5.645	3.854	3.774	-	1.867
S60S	158.2	117.4	59.95	53.96	34.94	23.89	12.73	7.944	7.752	-	3.226
N100S	327.6	236.6	113.2	100.9	62.51	41.09	20.51	12.19	11.86	-	4.538
S200S	639.5	451.3	204.7	180.8	107.9	68.64	32.38	18.44	17.91	-	6.350
S600S	2141	1454	599.5	521.2	291.0	173.9	73.68	36.65	37.36	-	11.41
S2000S	7929	5214	1982	1699	889.6	499.7	188.8	90.27	86.83	-	22.16
S8000S	34390	22088	7883	6684	3319	1767	598.1	258.7	247.8	-	51.24
N30000S	128417	80886	-	-	-	5595	-	-	-	-	111.6
N600	-	1454	-	-	-	173.9	-	-	-	15.40	-
N2000	-	5222	-	-	-	499.2	-	-	-	31.29	-

 VOTRE PARTENAIRE LABORATOIRE	<b>Agence Nord :</b> Plug N'Work - Campus Effiscience 2 Rue Jean Perrin - Bât D 14460 Colombelles Tél: 02.31.34.50.74 Fax: 02.31.34.55.17
	<b>Agence Sud :</b> Hôtel d'Entreprises de La Croix Rouge - Lot A4 10 Av de la Croix Rouge - 84000 Avignon Tél : 04.90.27.17.95 Fax : 04.90.27.17.52
<a href="http://www.deltalabo.fr">www.deltalabo.fr</a>	