

## Rotational Viscosity Standards

Product Code	Type	mPa.s (cP) 20°C	mPa.s (cP) 25°C	Pack Size
VIS-RT5-600	RT5	5.479	5.013	600mL
VIS-RT10-600	RT10	11.03	10.04	600mL
VIS-RT50-600	RT50	55.61	50.49	600mL
VIS-RT75-600	RT75	83.61	75.53	600mL
VIS-RT100-600	RT100	111.7	100.7	600mL
VIS-RT250-600	RT250	277.4	251.5	600mL
VIS-RT350-600	RT350	395.0	357.0	600mL
VIS-RT500-600	RT500	555.2	502.6	600mL
VIS-RT1K-600	RT1000	1107	1002	600mL
VIS-RT5K-600	RT5000	5523	4984	600mL
VIS-RT10K-600	RT10000	11265	10140	600mL
VIS-RT12K-600	RT12000	13718	12371	600mL
VIS-RT30K-600	RT30000	33291	30036	600mL
VIS-RT60K-600	RT60000	67075	60445	600mL
VIS-RT100K-600	RT100000	110129	99331	600mL

*All values stated are nominal - 28/08/2011*