

## Fuels Testing Standards

| Product Code | Test Method                          | Product Nature      | Certified Value | Pack Size |
|--------------|--------------------------------------|---------------------|-----------------|-----------|
| CRM-ADKR     | Acidity ASTM D 3242                  | Kerosene            | 0.0021 mg KOH/g | 250mL     |
| CRM-APKR     | Aniline Point ASTM D 611             | Jet Aviation fuel   | 62.30°C         | 250mL     |
| CRM-ACKR     | Aromatics Content ASTM D 5580        | Gasoline            | 28.25%          | 250mL     |
| CRM-BEGA     | Benzene Content ASTM D 4053 / D 5580 | Gasoline            | 0.827 % Vol.    | 250mL     |
| CRM-CNGO     | Cetane Number ASTM D 613             | Diesel              | 52.6            | 1000mL    |
| CRM-FIKR     | FIA Aromatics ASTM D 1319            | Kerosene            | 15.18 % vol     | 250mL     |
| CRM-OMGA     | Motor Octane Number ASTM D 2700      | Gasoline (Unleaded) | 86.9            | 1000mL    |
| CRM-ORGA     | Res. Octane Number ASTM D 2699       | Gasoline (Unleaded) | 98.0            | 1000mL    |
| CRM-VPGA     | Reid Vapour Pres. ASTM D 4953A, 323A | Gasoline            | 82.6 kPa        | 250mL     |

*All values stated are nominal - 12/07/2011*

### **Agence Nord:**

ZA Object'Ifs Sud - Lot A3  
6 Allée Emilie du Châtelet  
14123 Ifs  
tél : 02.31.34.50.74  
fax : 02.31.34.55.17



### **Agence Est:**

Parc Club des Tanneries  
2 Rue de la Faisanderie  
67380 Lingolsheim  
tél : 03.88.04.01.81  
fax : 03.68.93.01.52

### **Agence Sud:**

Bât Le Venango. 392 Rue Jean Dausset  
AGROPARC - BP11575  
84916 Avignon Cédex 9  
tél : 04.90.27.17.95 fax : 04.90.27.17.52

**www.deltalabo.fr**  
**info@deltalabo.fr**