

Reusable Threaded Centrifuge Tube, Graduated with White Scale, Unattached Cap

Large-capacity, 100 mL tube with 28-410 GPI screw thread closure. Cap supplied has PTFE-faced liner. Graduated and calibrated "to contain". Scale and marking spot are durable white ceramic enamel. Dimensions from 8" oil and weathering tube used in ASTM D2158.

Cap. (mL)	Subdivision (mL)	Tolerance ±(mL)	O.D. x Approx. Length (mm)	Shelf-Pack Qty.	Case Qty.
100	0-1 in 0.05	0-0.1 - ±0.02 Above 0.1-0.30 - 0.03 Above 0.3-1 - 0.05	37 x 203	6	12
	1-3 in 0.1	Above 1-3 - 0.10			
	3-6 in 0.2	Above 3-5 - 0.20			
	6-10 in 0.5	Above 5-10 - 0.50			
	10-100 in 1	Above 10 - 1.00			

Agence Nord:

ZA Object'Ifs Sud - Lot A3
6 Allée Emilie du Châtelet
14123 Ifs
tél : 02.31.34.50.74
fax : 02.31.34.55.17

Agence Est:

Parc Club des Tanneries
2 Rue de la Faisanderie
67380 Lingolsheim
tél : 03.88.04.01.81
fax : 03.68.93.01.52



Agence Sud:

Bât Le Venango. 392 Rue Jean Dausset
AGROPARC - BP11575
84916 Avignon Cédex 9
tél : 04.90.27.17.95 fax : 04.90.27.17.52

www.deltalabo.fr
info@deltalabo.fr



