

Robinson Centrifuge Model 872AT4



Features:

- Solid cast aluminum housing designed to withstand harsh field conditions
- Variable speed motor control allows for slow, smooth start.
- Modular component design allows for easy removal and replacement of motor, rheostat, switches, and heaters without unmounting the centrifuge.
- Pre-heaters on the front of the machine are easily accessible to the operator
- High-efficiency bowl heaters maintain sample temperature for the duration of testing.
- 15-minute analog timer automatically stops motor when testing is complete.
- Spring loaded lid lock keeps lid closed during operation.
- Center hole in lid gives safe and easy access to check RPM's with a handheld tachometer during operation.

Power Source	230 VAC, 50/60 Hz
Bowl Heater	Yes
Pre-Heater	4-Place
Tube Style	100mL Short Cone
Trunnion Head	4-Place
Timer	15-Minute Analog
Tachometer	N/A
Temperature Controller	N/A
Class 1, Division 2	No
Maximum RPM (RCF)	1700±100 (675±82)
Product Dimensions (in)	23.5L x 20W x 11.5H
Product Weight (lbs)	72



Agence Nord:
ZA Object'ifs Sud - Lot A3
6 Allée Emilie du Châtelet
14123 Ifs
tél : 02.31.34.50.74
fax : 02.31.34.55.17

Agence Est:
Parc Club des Tanneries
2 Rue de la Faisanderie
67380 Lingolsheim
tél : 03.88.04.01.81
fax : 03.68.93.01.52



Agence Sud:
Bât Le Venango, 392 Rue Jean Dausset
AGROPARC - BP11575
84916 Avignon Cédex 9
tél : 04.90.27.17.95 fax : 04.90.27.17.52

www.deltalabo.fr
info@deltalabo.fr