

Agence Nord:
ZA Object'ifs Sud - Lot A3
6 Allée Emilie du Châtelet
14123 Ifs
tél : 02.31.34.50.74
fax : 02.31.34.55.17



Agence Est:
Parc Club des Tanneries
2 Rue de la Faisanderie
67380 Lingolsheim
tél : 03.88.04.01.81
fax : 03.68.93.01.52

Agence Sud:
Bât Le Venango, 392 Rue Jean Dausset
AGROPARC - BP11575
84916 Avignon Cédex 9
tél : 04.90.27.17.95 fax : 04.90.27.17.52

www.deltalabo.fr
info@deltalabo.fr

TANAKA
Petroleum Testing & Beyond

EDXRF Sulfur Analyzer

rx-460SH

NEW

Compliant to ASTM D4294-03

Limit of Detection 15 mg/kg

Portable and easy to use



Specifications

Test Method	ASTM D4294-03, ISO 8754, JIS K 2541-4, etc.
Measuring Range	15 mg/kg (LOD) to 6.00 wt%
Size (W x D x H)	430 mm x 410 mm x 180 mm
Weight	16.2 kg
Power Consumption	AC 100 V to 240 V 100 VA (max) 50 / 60 Hz

We reserve the right of changes without prior notice.



202509P

Tanaka Scientific Limited

7-10-3, Ayase, Adachi-ku, Tokyo, 120-0005 JAPAN
TEL : +81-3-3620-1711 E-mail: overseas-group@tanaka-sci.com
FAX: +81-3-3620-1713 URL : www.tanaka-sci.com