

**Agence Nord:**

ZA Object'Ifs Sud - Lot A3  
 6 Allée Emilie du Châtelet  
 14123 Ifs  
 tél : 02.31.34.50.74  
 fax : 02.31.34.55.17

**Agence Est:**

Parc Club des Tanneries  
 2 Rue de la Faisanderie  
 67380 Lingolsheim  
 tél : 03.88.04.01.81  
 fax : 03.68.93.01.52

**Agence Sud:**

Bât Le Venango. 392 Rue Jean Dausset  
 AGROPARC - BP11575  
 84916 Avignon Cédex 9  
 tél : 04.90.27.17.95 fax : 04.90.27.17.52

[www.deltalabo.fr](http://www.deltalabo.fr)  
[info@deltalabo.fr](mailto:info@deltalabo.fr)



## Osmolality standards

Osmolality Standards for calibration and control, compatible with all types of osmometers.

*Reagecon*

Extensive range: 50 to 3000 mOsmol/kg H<sub>2</sub>O, ISO 17025 certified. Values from 100 to 700 mOsmol/kg H<sub>2</sub>O according to European Pharmacopoeia 2.2.35 and US Pharmacopoeia 785. 12 ampoules of 5 ml each.

Type	Description	PK	Cat. No.
RE-OSM-50	50 mOsmol/kg H <sub>2</sub> O	12	<b>6.286 326</b>
RE-OSM-100	100 mOsmol/kg H <sub>2</sub> O	12	<b>6.286 324</b>
RE-OSM-200	200 mOsmol/kg H <sub>2</sub> O	12	<b>6.286 325</b>
RE-OSM-290	290 mOsmol/kg H <sub>2</sub> O	12	<b>6.280 786</b>
RE-OSM-300	300 mOsmol/kg H <sub>2</sub> O	12	<b>6.284 310</b>
RE-OSM-400	400 mOsmol/kg H <sub>2</sub> O	12	<b>7.970 554</b>
RE-OSM-500	500 mOsmol/kg H <sub>2</sub> O	12	<b>6.270 201</b>
RE-OSM-600	600 mOsmol/kg H <sub>2</sub> O	12	<b>6.286 327</b>
RE-OSM-700	700 mOsmol/kg H <sub>2</sub> O	12	<b>6.286 328</b>
RE-OSM-850	850 mOsmol/kg H <sub>2</sub> O	12	<b>6.286 334</b>
RE-OSM-900	900 mOsmol/kg H <sub>2</sub> O	12	<b>6.286 335</b>
RE-OSM-1000	1000 mOsmol/kg H <sub>2</sub> O	12	<b>6.286 329</b>
RE-OSM-1500	1500 mOsmol/kg H <sub>2</sub> O	12	<b>6.286 330</b>
RE-OSM-2000	2000 mOsmol/kg H <sub>2</sub> O	12	<b>6.286 331</b>
RE-OSM-2500	2500 mOsmol/kg H <sub>2</sub> O	12	<b>6.286 332</b>
RE-OSM-3000	3000 mOsmol/kg H <sub>2</sub> O	12	<b>6.286 333</b>