

Agence Nord:
ZA Object'Ifs Sud - Lot A3
6 Allée Emilie du Châtelet
14123 Ifs
tél : 02.31.34.50.74
fax : 02.31.34.55.17



Agence Est:
Parc Club des Tanneries
2 Rue de la Faisanderie
67380 Lingolsheim
tél : 03.88.04.01.81
fax : 03.68.93.01.52

Agence Sud:
Bât Le Venango, 392 Rue Jean Dausset
AGROPARC - BP11575
84916 Avignon Cédex 9
tél : 04.90.27.17.95 fax : 04.90.27.17.52

www.deltalabo.fr
info@deltalabo.fr

40-282-001 Automatic Sample Feeder (Robot for 6-station VST ECO)



Application

Automatic Sample Feeding as add-on to the VST ECO for fully automated VST testing.

Features

Sample Holder for up to 120 Samples. One set from 1-6 samples depending on test configuration are fed into the test chamber. After testing the old samples are removed. The feeding system is connected to the machine controller and is fully integrated in the measuring process.

With this robot system the Coesfeld VST ECO runs fully independent and safely over night or weekend. Including the robot system, the machine is not significantly larger than many oil based tester without automatisation.

Technical Data

Capacity 120 samples 10x10x4 mm

Dimension and Connection

Dimension (WxDxH)	1200 x 750 x 1150 mm (with VST ECO)
Weight	30 kg
Interface	Direct link to VST ECO (6-station)
Air	Direct link to VST ECO (6-station)