




61-400 Apparent Density Tester

Standards

DIN 53 466 - ISO 60 - ISO 171 – GB/T16362008 Method A



Application

Measuring system for determination the apparent density of free flow powder or granulate. The test can provide a gross measure of particle size and dispersion, which can affect material flow consistency and reflect packaging quantity.

Features

- Solid, warp resistant stand
- Height adjustable loading funnel
- Snap closure for loading funnel with recuperating spring
- Scraper made of stainless steel
- Calibrated measuring cup
- Drip cup for surplus material

Dimensions

| | |
|--------------------|---------------------------------------|
| Dimensions (WxDxH) | 260 x 160 x 210 mm |
| Weight | approx. 3.5 kg |
| Funnel | height: 115 mm, volume approx. 200 ml |
| Measuring cup | 100 ±0.5 ml |

| | | |
|---|--|---|
| <p>Agence Nord: ZA Objectifs Sud - Lot A3 6 Allée Emilie du Châtelet 14123 Ifs tél : 02.31.34.50.74 fax : 02.31.34.55.17</p> | | <p>Agence Est: Parc Club des Tanneries 2 Rue de la Faisanderie 67380 Lingolsheim tél : 03.88.04.01.81 fax : 03.68.93.01.52</p> |
| <p>Agence Sud: Bât Le Venango, 392 Rue Jean Dausset AGROPARC - BP11575 84916 Avignon Cédex 9 tél : 04.90.27.17.95 fax : 04.90.27.17.52</p> | | <p>www.deltalabo.fr info@deltalabo.fr</p> |