

**Agence Nord:**  
ZA Object'Ifs Sud - Lot A3  
6 Allée Emilie du Châtelet  
14123 Ifs  
tél : 02.31.34.50.74  
fax : 02.31.34.55.17




**Agence Est:**  
Parc Club des Tanneries  
2 Rue de la Faisanderie  
67380 Lingolsheim  
tél : 03.88.04.01.81  
fax : 03.68.93.01.52

**Agence Sud:**  
Bât Le Venango, 392 Rue Jean Dausset  
AGROPARC - BP11575  
84916 Avignon Cédex 9  
tél : 04.90.27.17.95 fax : 04.90.27.17.52

[www.deltalabo.fr](http://www.deltalabo.fr)  
[info@deltalabo.fr](mailto:info@deltalabo.fr)

## 50-014 / 50-014.1 – Kofler Heating Bench



### Application

- Determination of melting point
- Determination of softening point of synthetic resins
- Determining the degree of purity of substances
- Presence and behaviour of water of crystallisation
- Observing decomposition processes and volatility
- Determining eutectic temperatures
- Determining the softening or fusion point of synthetic resin powders

### Features

Incl. test- and gauge substance set (8 pieces)

### Technical Data

<b>Mains:</b>	230 V / 50/60 Hz (item no. 50-014) 115 V / 50/60 Hz (item no. 50-014.1)
<b>Power:</b>	max. 100 W
<b>Temperature range:</b>	+50...+260°C
<b>Temperature scale graduations:</b>	2°C
<b>Measuring accuracy in the calibrated range:</b>	±1 °C
<b>Dimensions:</b>	approx. 430 x 100 x 140 mm
<b>Weight:</b>	approx. 2.5 kg